> podo+ech® gels

gel protection and insole range





A. Algeo Ltd Sheridan House Bridge Industrial Estate Speke Hall Road Liverpool, L24 9HB, UK

| loe separators and loe Caps | rage 3 |
|---|---------|
| A selection of gel separators and caps designed to relieve toe pressure and provide comfort M-Gel Toe Separators • M-Gel Toe Spreaders • M-Gel Universal Toe Crest • M-Gel Toe/Finger Caps | |
| Silicone Digits | Page 4 |
| Interdigital padded products for cushioning and pain relief M-Gel Digital Pads • M-Gel Corn Pads • M-Gel Digital Cap • Digital Pad Strips • Coated Digital Tubes | |
| Metatarsal Cushioning | Page 5 |
| Providing protection for the metatarsal area S-Gel Metatarsal Pads • Forefoot Relief Pads • M-Gel Metatarsal Straps • M-Gel Forefoot Protectors | |
| Bunion Guards | Page 5 |
| A selection of products designed specifically for the protection against bunions M-Gel Thin Dress Bunion Reliefs • M-Gel Big Toe Bunion Guards • M-Gel 5th Toe Bunion Guards | |
| Ankle Protection | Page 6 |
| Providing protective cushioning and pressure relief around the ankle Achilles Protection Sleeves • Ankle Protection Sleeves • Dorsum Protection Sleeves | |
| Gel Squares and Silicone Insoles | Page 7 |
| Designed to relieve pressure and absorb shock M-Gel Squares • M-Gel Dots Sheets | |
| Standard Dress Insoles | Page 7 |
| A range of dress insoles designed to provide relief in the tightest of shoes S-Gel Dress Insoles S-Gel 3/4 Length Insoles S-Gel Thin Insoles S-Gel Thin Insoles Met Raise | |
| Heel Cups | Page 10 |
| A selection of heel cups designed to provide shock absorption and cushioning for the heel Heel Cups with Removable Spurs • Heel Cup with Shock Pattern • Dress Heels • Dual Density Heel Cups | |
| Heel Wedges | Page 1 |
| Comfort and pain relief designed for all types of foot S-Gel Posted Heel Wedges • S-Gel Heel Wedge with Cover • S-Gel Medium Density Heel Wedge | |
| Moisturising Accessories | Page 12 |
| A selection of socks to moisturise and rejuvenate dry and tired feet M-Gel Moisturising Rooties & M-Gel Moisturising Heel Sleeve & M-Gel Moisturising Vented Heel Sleeve | |

About the range

Designed to absorb shock as well as relieve and redistribute pressure, the Podotech gel range offers a large selection of gel products from toe separators and pads to insoles. The range is manufactured from high performance medical grade gel designed to provide optimal functionality and best value.

You can count on the Podotech range of gel products to give consistent quality, which means that reliable performance from every product, every time.

Podotech insoles stay in place for unsurpassed comfort and gel coatings stand up to repeated wear and washings. There are two types of gel in the range which are designed to meet both performance and cost parameters.



Washable and reusable - can stand up to repeated wash



Highly Durable - will not flatten or lose its shape like foam products



M-Gel releases moisture and vitamir on contact with skir



S.C.o.

Silicone Gel is universally recognised and accepted by the medical community as an excellent solution to podiatric cushioning. Silicone is considered to be the traditional top-of-the-line material for contact with the skin. Silicone gel is extremely durable and hypoallergenic.

M-Gel

Mineral Oil Gel is a thermoplastic elastomer that is 100% medical grade and non-allergenic. The gel is an enriched formula that releases moisture and vitamins on contact, softening and soothing the driest of skin.





Podotech Gel products are manufactured using Polygel. FabricSmart is the trademark of Gel Concepts LLC

M-Gel



M-Gel Toe Separator

| Code | Size |
|--------|-------------|
| PT0001 | Small |
| PT0002 | Large |
| PT0003 | Extra Large |

- Traditional anatomical shape designed for standard footwear to minimise toe pressure
- Fits comfortably between toes to absorb pressure
- Helps to reduce friction between toes with vitamin enriched mineral oil gel
- Relieves pain from soft corns, overlapping toes, nail problems and bunions



M-Gel Toe Spreader

| Code | Size |
|--------|--------|
| PT0012 | Small |
| PT0015 | Medium |
| PT0018 | Large |

- For straightening and aligning toes and relieving pressure form bunions
- Vitamin enriched gel cushions and protects the skin



M-Gel Universal Toe Crest

| Code | Size |
|--------|--------|
| PT0043 | Medium |

- To ease discomfort for hammer, claw, and mallet toe
- Comfortable gel exudes medical grade mineral oil to soothe the skin
- Flexible toe loop easily fits over toe to stay in place



M-Gel Toe Cap

| Code | Size |
|--------|---------------|
| PT0071 | Small/Medium |
| PT0073 | Large/X Large |
| PT0075 | XX Large |
| | |

- Traditional shape for relieving pressure and friction on toes or fingers
- Superior cushioning and relief for corns, calluses, hammer toes and partially amputated digits
- Comfort vent for breathability
- Clear material allows discrete use in open sandals

silicone digits



M-Gel Digital Pads

| Code | Size |
|--------|---------------|
| PT0076 | Mini |
| PT0077 | Small/Medium |
| PT0078 | Large/X Large |

- Thicker, round shaped gel padding for increased cushioning and comfort
- A soft, comfortable and stretchable slipon device
- Relieves pain from corns, blisters, hammer toes, calluses and digital trauma
- Gel coating softens, soothes and protects the skin

M-Gel



M-Gel Corn Pads

| Code | Size |
|--------|---------------|
| PT0082 | Small/Medium |
| PT0083 | Large/X Large |

- Traditional product for protecting toes or fingers
- Flexible fabric coated one half with soft gel
- Gel design for maximum coverage
- Relieves pain from corns, blisters, hammer toes, calluses and digital injuries
- Gel coating softens, soothes and protects the skin





M-Gel Digital Caps

| Code | Size |
|--------|---------------|
| PT0087 | Small/Medium |
| PT0088 | Large/X Large |
| PT0089 | XX Large |
| | |

- Unique seamless construction provides comfort and optimal toe/finger protection
- Gel softens, soothes and protects
- Superior cushioning, protection and relief for corns, calluses, nails and cuticles, hammer toes and partially amputated digits
- Improves skin elasticity, as well appearance of old and new scars



M-Gel Digital Pad Strips

| Code | Size |
|--------|---------------|
| PT0097 | Mini |
| PT0098 | Small/Medium |
| PT0099 | Large/X Large |
| | |

- Protection for toes/fingers from pressure, shock and shear
- Uniform round dots of soft, comfortable gel on each strip.
- Can be easily cut to place padding where needed
- Relieves pain from corns, calluses, nail problems, dry or damaged cuticles



M-Gel Coated Digital Tubes

| Code | Size |
|--------|--------|
| PT0104 | Mini |
| PT0105 | Narrow |
| PT0106 | Wide |

- Cut to size as and when you wish
- Fits comfortably between toes to absorb pressure
- Reduces friction and irritation between toes with vitamin enriched gel
- Money saver for heavy users of digital corn pads



M-Gel Thin Forefoot Pad

| Code | Size |
|--------|---------------|
| PT0151 | Small/Medium |
| PT0152 | Large/X Large |
| | |

- Ideal for the protection of metatarsal heads
- Helps to relieve burning sensation under the metatarsal heads
- Soft, gel pad with flexible, thin elastic loop attachable over the second or third toe
- Fabric cover prevents sliding in footwear

metatarsal cushioning



S-Gel Metatarsal Pad & Toe Separator

| Code | Size |
|--------|-------------|
| PT0125 | Small Right |
| PT0126 | Small Left |
| PT0127 | Large Right |
| PT0128 | Large Left |

S-Gel Metatarsal Pad with Loop

| Code | Туре |
|--------|-------|
| PT0129 | Left |
| PT0130 | Right |

- Soft gel cushions while effectively absorbing shock and shear forces
- Protects the metatarsal heads
- Ideal for fat pad atrophy and metatarsalgia



Thermoplastic Forefoot Relief Pad with Cover

| Code | Size |
|--------|----------|
| PT0134 | One Size |

- Cushions the ball of the foot to help prevent aches or pains
- Designed with a shock absorbing pattern which holds the product in place
- Recommended for high heels, sandals or boots
- Durable to provide all day cushioning
- Anti-microbial top cover absorbs odour and prevents bacterial growth



Thermoplastic Forefoot Pad

| Code | Size |
|--------|----------|
| PT0133 | One Size |

- Cushions the ball of the foot to help prevent aches or pains
- Designed with a shock absorbing pattern which holds the product in place
- Recommended for high heels, sandals or boots
- Durable to provide all day cushioning
- Transparent to match any shoe



M-Gel Metatarsal Strap

| Code | Size |
|--------|---------------------|
| PT0137 | Small/Medium Right |
| PT0138 | Small/Medium Left |
| PT0139 | Large/X Large Right |
| PT0140 | Large/X Large Left |

- For protection of the metatarsal heads from pressure and shear forces
- Helps to relieve burning sensation under the metatarsal heads and other pain from calluses, Morton's Neuroma and fat pad atrophy
- Flexible, durable outer fabric for comfort and conformability
- Redistributes and eases impact pressure from walking, running or carrying excess weight



M-Gel Forefoot Protector

| Code | Size |
|--------|-------|
| PT0153 | Small |
| PT0154 | Large |

- Constructed from a thin layer of gel covered on both sides with fabric
- Provides complete toe protection
- Extends to beneath the metatarsal heads
- Worn comfortably in every day footwear and under hosiery



M-Gel Thin Dress Bunion Relief

| Code | Size |
|--------|-------|
| PT0168 | Small |
| PT0169 | Large |

- 3mm soft flexible footpad for protecting bunions by absorbing shock and shear forces
- soft mineral oil gel soften and comforts painful and sensitive bunions on direct contact with the
- Secures comfortably over the big toe and around the foot to prevent sliding

ankle protection



M-Gel Big Toe Bunion Guard

| Code | Size |
|--------|----------|
| PT0176 | One Size |

- Ideal for protection of a hallux bunion from friction or pressure
- Soft flexible shell cushions and comforts tender and sensitive skin and surrounding tissue
- Comfortable toe loop slides over the big toe to hold the product in place
- Thin profile design



M-Gel 5th Toe Bunion Guard

| Code | Size |
|--------|----------|
| PT0177 | One Size |

- Ideal for protection of the 5th digit
- Soft flexible shell cushions and comforts tender and sensitive skin and surrounding tissue
- Comfortable toe loop slides over the pinky toe to hold the product in place
- Thin profile design



M-Gel Achilles Protection Sleeve

| Code | Size |
|--------|-------|
| PT0178 | Small |
| PT0179 | Large |

- Ideal for relief from pump bones (Haglund's deformity) or sore, fatigued or injured achilles tendon
- Soft gel pad protects the posterior of the heel from friction whilst gently compressing and reducing tension around the heel area
- Designed for slim profile in footwear



M-Gel Ankle Protection Sleeve

| Code | Size |
|--------|----------|
| PT0180 | One Size |

- For protective cushioning around ankle bones and relief of pressure building up during high impact sports
- Soft, comfortable pad provides gentle compression around ankle bones whilst reducing abrasion and shear forces over the
- Ideal ankle protection in boots and braces



M-Gel Dorsum Protection Sleeve

| Code | Size |
|--------|-------|
| PT0184 | Small |
| PT0185 | Large |

- For protection from the top of the foot to the shin by reducing pressure, shear forces and friction
- Popular product used as an interface between the foot and the boot
- Comfortable
 elastic sleeve holds
 moulded,
 anatomically
 designed gel
 shape in place on
 top of the foot



M-Gel Achilles & Dorsum Protection Sleeve

| Code | Size |
|--------|-------|
| PT0186 | Small |
| PT0187 | Large |
| | |

- Protection for both the top of the foot and Achilles tendon by reducing pressure, shear forces and friction from boots
- Protects from abrasions caused by athletic footwear without added
- Comfortable gel sleeve holds moulded, anatomically designed gel shape in place on top of the foot and supporting the Achilles tendon

gel squares & insoles



M-Gel Square

| Code | Size | Туре |
|--------|-------------|------------------|
| PT0194 | 101 x 101mm | Without Adhesive |
| PT0195 | 101 x 101mm | With Adhesive |

- Soft formulation for anywhere on the body to offload pressure and absorb shock or shear forces
- 6mm thick pure mineral oil gel available with or without pressure sensitive adhesive backing
- May help improve the appearance of old or new scars



M-Gel Dots Sheet

| Code | Size | Туре |
|--------|------|---------------|
| PT0196 | 25mm | With Adhesive |

- Soft formulation for anywhere on the body to offload pressure and absorb shock or shear forces
- Peel-and-stick adhesive backing for secure placement anywhere
- Mineral oil moisturises skin at point of contact
- May help improve the appearance of old or new scars



S-Gel Thin Dress Insole with Metatarsal Raise

| Code | Size | |
|--------|--------|--|
| PT0208 | Small | |
| PT0209 | Medium | |
| PT0210 | Large | |
| | | |

- For absorbing shock, enhanced comfort and forefoot pressure relief while wearing high-heeled or tightfitting shoes
- Thin design fits into most women's dress footwear
- Soft metatarsal raise redistributes pressure on forefoot for relief and and comfort while providing added support and stability



S-Gel Thin Dress Insole with Metatarsal Raise - Covered

| | Code | Size | |
|--|--------|--------|--|
| | PT0211 | Small | |
| | PT0212 | Medium | |
| | PT0213 | Large | |
| | | | |

fabricsmart

- For absorbing shock, enhanced comfort and forefoot pressure whilst wearing high heeled and tight fitting shoes
- Thin design fits into most women's dress footwear
- Anti-microbial cover

About anti-microbial coverings

Anti-microbial top covers feature Eco-Fibre Shield technology that is environmentally and user friendly. No metals means no worries.

Advanced Fabric Smart technology greatly increases the practical value and comfort of silicone orthotic shoe inserts. All of Fabric Smart antimicrobial fabrics provides the following benefits:

- Antibacterial, Antifungal and Deodourising Our proprietary treatment is infused into the fabric so it will not rub off or migrate onto the skin
- Breathable and dry Absorbs and wicks away moisture and perspiration fast, while providing excellent ventilation, keeping the skin cool and dry



- 4-Way stretch fabric Will not deform or flatten after long usage
- **Durable** Fabric will not undergo deformation or alter thickness after long usage
- Washable Maintains effectiveness after washing

FabricSmart is the trademark of Gel Concepts LLC

| shoe size conversion guide | | | | | | | | | | | | |
|----------------------------|------|------|------|----|----|----|----|----|------|----|----|----|
| UK Men | 3 | 31/2 | 41/2 | 5½ | 6 | 7 | 8 | 9 | 91/2 | 10 | 11 | 12 |
| UK Women | 31/2 | 4 | 5 | 6 | 6½ | 7½ | 8 | 9 | | | | |
| EU | 36 | 37 | 38 | 39 | 40 | 41 | 42 | 43 | 44 | 45 | 46 | 47 |

dress insoles



S-Gel Dress Insole With Wide Metatarsal Pad - Covered

| Code | Size | |
|--------|--------|--|
| PT0219 | Small | |
| PT0220 | Medium | |
| PT0221 | Large | |

- For absorbing shock, enhanced comfort and forefoot pressure whilst wearing highheeled and tightfitting shoes
- Thin design fits into most women's and men's dress footwear
- Soft, wide metatarsal pad provides maximum comfort while easing pressure on the forefoot
- Anti-microbial cover

fabricsmart



S-Gel Dress Insole With Wide Metatarsal Pad

| Code | Size | |
|--------|--------|--|
| PT0214 | Small | |
| PT0215 | Medium | |
| PT0216 | Large | |

- For absorbing shock, enhanced comfort and forefoot pressure relief while wearing high-heeled or tight-fitting shoes
- Thin design fits into most women's and men's dress footwear
- Soft, wide metatarsal pad provides maximum comfort while easing pressure on the forefoot



S-Gel Thin Dress Insole - Covered

| Code | Size |
|--------|--------|
| PT0203 | Small |
| PT0204 | Medium |
| PT0205 | Large |

- Thin design fits into most women's dress footwear
- Soft metatarsal raise redistributes pressure on forefoot for relief and and comfort while providing added support and stability
- Anti-microbial cover





S-Gel Thin Insole With Reliefs

| Code | Size | EU | UK (Men's) |
|--------|---------|----|------------|
| PT0232 | Small | 35 | 3 |
| PT0233 | Medium | 38 | 41/2 |
| PT0234 | Large | 41 | 7 |
| PT0235 | X Large | 44 | 9 |
| | | | |

- Full foot coverage with dual durometer construction for pressure point relief and comfort
- Provides extra shock absorption under the heel and forefoot, with a soft metatarsal pad to minimise pressure on metatarsal heads
- Silicone surface features vent holes for enhanced circulation



S-Gel Thin Insole With Reliefs - Covered

| Code | Size | EU Size | UK (Men's) |
|--------|---------|---------|------------|
| PT0236 | Small | 35 | 3 |
| PT0237 | Medium | 38 | 41/2 |
| PT0238 | Large | 41 | 7 |
| PT0239 | X Large | 44 | 9 |

- Full foot coverage with dual durometer construction for pressure point relief and comfort
- Provides extra shock absorption under the heel and forefoot, with soft a metatarsal pad to minimise pressure
- Silicone surface features vent holes for enhanced circulation
- Anti-microbial cover





S-Gel Flat Insole - Covered

| Size | EU Size | UK (Men's) |
|--------|-----------------|-----------------------|
| Small | 38 | 41/2 |
| Medium | 43 | 9 |
| Large | 46 | 11 |
| | Small Medium | Small 38 Medium 43 |

- Long lasting basic shape fits on top of built in insoles for complete gel comfort from toe box to heel
- Absorbs shock and shear forces
- Fits into women's and men's dress shoes
- Anti-microbial cover



silicone insoles



S-Gel Flat Insole with Metatarsal Raise and Toe Crest - Covered

| | Code | Size | EU Size | UK (Men's) |
|--|--------|---------|---------|------------|
| | PT0246 | Small | 34 | 3 |
| | PT0247 | Medium | 37 | 31/2 |
| | PT0248 | Large | 40 | 6 |
| | PT0249 | X Large | 43 | 9 |

- Corrective product for relieving heel pain and improving stability
- Wide, comfortable metatarsal raise
- Toe crest helps straighten hammer toes and ease discomfort from other toe problems
- Low profile design made for a comfortable fit in today's footwear
- Anti-microbial cover

fabricsmart



S-Gel Insole With Soft Reliefs

| Code | Size | EU Size | UK (Men's) |
|--------|-----------|---------|------------|
| PT0250 | Small | 36 | 3 |
| PT0251 | Medium | 38 | 41/2 |
| PT0252 | Large | 40 | 6 |
| PT0253 | X Large | 42 | 8 |
| PT0387 | XX Large | 44 | 9 |
| PT0388 | XXX Large | 46 | 11 |

- Product for relieving heel pain and improving stability while redistributing pressure away from common stress points
- Wide, comfortable metatarsal raise
- Thicker design works well in active wear shoes or boots
- Silicone surface features vent holes for enhanced circulation



S-Gel Insole With Soft Reliefs And Metatarsal Raise

| Code | Size | EU Size | UK (Men's) |
|--------|-----------|---------|------------|
| PT0260 | Small | 36 | 3 |
| PT0261 | Medium | 38 | 41/2 |
| PT0262 | Large | 40 | 6 |
| PT0263 | X Large | 42 | 8 |
| PT0264 | XX Large | 44 | 9 |
| PT0265 | XXX Large | 46 | 11 |
| | | | |

- Dual durometer construction for extra protection of heel, mid-foot and forefoot areas
- Provides extra shock absorption under the metatarsal heads with a soft metatarsal pad to minimise pressure
- Silicone surface features vent holes for enhanced circulation



S-Gel Insole With Soft Reliefs - Covered

| Code | Size | EU Size | UK (Men's) |
|--------|-----------|---------|------------|
| PT0254 | Small | 36 | 3 |
| PT0255 | Medium | 38 | 41/2 |
| PT0256 | Large | 40 | 6 |
| PT0257 | X Large | 42 | 8 |
| PT0258 | XX Large | 44 | 9 |
| PT0259 | XXX Large | 46 | 11 |
| | | | |

- Product for relieving heel pain and improving stability while redistributing pressure away from common stress points
- Wide, comfortable metatarsal raise
- Thicker design works well in active wear shoes or boots
- Silicone surface features vent holes for enhanced circulation
- Anti-microbial cover

fabricsmart



S-Gel Insole With Soft Reliefs, Metatarsal Raise - Covered

| Code | Size | EU Size | UK (Men's) |
|--------|-----------|---------|------------|
| PT0267 | Small | 36 | 3 |
| PT0268 | Medium | 38 | 41/2 |
| PT0269 | Large | 40 | 6 |
| PT0270 | X Large | 42 | 8 |
| PT0271 | XX Large | 44 | 9 |
| PT0389 | XXX Large | 46 | 11 |
| | | | |

fabricsmart

- Dual durometer construction for extra protection of heel, mid-foot and forefoot areas
- Provides extra shock absorption under the metatarsal heads with a soft metatarsal pad to minimise pressure
- Silicone surface features vent holes for enhanced circulation
- Anti-microbial cover



S-Gel Deep 3/4 Length Insole

| Code | Size | EU Size | UK (Men's) |
|--------|----------|---------|------------|
| PT0222 | Small | 36 | 3 |
| PT0223 | Medium | 38 | 41/2 |
| PT0224 | Large | 40 | 6 |
| PT0225 | X Large | 42 | 8 |
| PT0226 | XX Large | 44 | 9 |

- '1 Spot' 3/4 insoles also available on request.
- Dual durometer construction for extra protection of rear and mid-foot tissues without crowding the shoe
- Designed for flat shoes and trainers
- High medial and lateral walls cradle soft rear foot tissues
- Raised wide pad of softer silicone in the metatarsal area relieves pressure

heel cups



S-Gel 3/4 Length Insole - Covered

| Code | Size | EU Size | UK (Men's) |
|--------|----------|---------|------------|
| PT0227 | Small | 36 | 3 |
| PT0228 | Medium | 38 | 41/2 |
| PT0229 | Large | 40 | 6 |
| PT0230 | X Large | 42 | 8 |
| PT0231 | XX Large | 44 | 9 |

- Dual durometer construction for extra protection of rear and mid-foot tissues without crowding the shoe
- High medial and lateral walls cradle soft rear foot tissues
- Raised wide pad of softer silicone in the metatarsal area relieves pressure in the metatarsal heads
- Anti-microbial cover





S Gel

S-Gel Heel Pad with Removable Spur

| Code | Size |
|--------|--------|
| PT0283 | Small |
| PT0284 | Medium |
| PT0285 | Large |
| | |

- Adjustable with spur plug for shock absorption or removed for spur relief
- Medium density silicone efficiently absorbs impact and recovers quickly
- Improved smaller spur pad fits correctly under the heel aperture



S-Gel Heel Pad with Removable Spur - Covered

| Code | Size |
|--------|--------|
| PT0288 | Small |
| PT0289 | Medium |
| PT0290 | Large |

- Adjustable with spur plug for shock absorption or removed for spur relief
- Medium density silicone efficiently absorbs impact and recovers quickly
- Improved smaller spur pad fits correctly under the heel aperture
- Will not deform or degrade with use
- Anti-microbial cover





S-Gel Heel Cup With Shock Pattern - Covered

| Code | Size |
|--------|--------|
| PT0296 | Small |
| PT0297 | Medium |
| PT0298 | Large |

- Protection for general use and high-impact activities
- Waffle pattern in contact with the body maximises shock absorption
- Thin heel fits comfortably in most footwear
- Will not deform or degrade with use
- Anti-microbial cover





S-Gel Dress Heel With Spur Relief

| Code | Size |
|--------|--------|
| PT0299 | Small |
| PT0300 | Medium |
| PT0301 | Large |
| | |

- Spur relief in shallow fitting dress footwear
- Thin, flexible silicone heel wraps around the inside of the shoe to cradle the heel
- Design helps to secure the heel in place
- Works where no other heel support will fit



S-Gel Heel Cup

| Code | Size |
|--------|---------|
| PT0304 | Small |
| PT0305 | Medium |
| PT0306 | Large |
| PT0310 | X Large |
| | |

- Dual durometer for maximum comfort while standing, walking or running
- Medium durometer silicone in main body efficiently absorbs impact and recovers quickly
- Softer silicone blue zone under heel centre maximises shock absorption
- Will not deform or degrade with use

heel wedges



S-Gel Heel Cup - Covered

| Code | Size |
|--------|--------|
| PT0307 | Small |
| PT0308 | Medium |
| PT0309 | Large |

- Medium dual durometer for maximum comfort while standing, walking or running
- Softer silicone blue zone under heel centre maximises shock absorption and comfort at point of impact
- Will not deform or degrade with use
- Anti-microbial cover

fabricsmart Eco-Fiber Shield



S-Gel Spur Cup - Medial Side

| Code | Size |
|--------|--------|
| PT0315 | Small |
| PT0316 | Medium |
| PT0317 | Large |

- Medium dual durometer for medial-side or spur/fascia relief and comfort while standing, walking or running
- Softer silicone blue zone under tender areas provides increased comfort and pain relief
- Retains shape will not deform or degrade with use



S-Gel Spur Cup Medial Side - Covered

| Code | Size |
|--------|--------|
| PT0320 | Small |
| PT0321 | Medium |
| PT0322 | Large |



- Dual durometer cup for medial-side or spur/fascia relief and comfort
- Softer silicone blue zone under tender areas provides increased comfort and pain relief
- Will not deform or degrade with use
- Anti-microbial cover



S-Gel

S-Gel Dual Density Posted Heel Wedge

| Code | Size |
|--------|--------|
| PT0331 | Small |
| PT0332 | Medium |
| PT0333 | Large |
| | |

- Dual durometer wedge for medialside or spur/fascia relief for users with flat feet, high arches or knee conditions
- Softer silicone medial blue zone increases relief from heel spur conditions
- Posted design minimises pronation, flat feet or general valgus knee conditions
- Will not deform or degrade with use



S Gel

S-Gel Posted Heel Wedge - Covered

| Code | Size |
|--------|--------|
| Code | Size |
| PT0334 | Small |
| PT0335 | Medium |
| PT0336 | Large |
| | |

- Dual durometer for medial-side or spur/fascia relief
- Softer silicone increases comfort and pain relief from heel spur conditions
- Posted design minimises pronation, flat feet or general valgus knee conditions
- Anti-microbial cover

fabricsmart



S-Gel Heel Wedge

| Code | Size |
|--------|--------|
| PT0337 | Small |
| PT0338 | Medium |
| PT0339 | Large |

- Medium dual durometer for maximum comfort while standing, walking or running
- Softer silicone blue dot under heel centre maximises shock absorption and comfort at point of impact
- Will not deform or degrade with use

moisturising accessories



S-Gel Heel Wedge - Covered

| Code | Size |
|--------|--------|
| PT0340 | Small |
| PT0341 | Medium |
| PT0342 | Large |
| | |

- Medium dual durometer for maximum comfort while standing, walking or running
- Softer silicone blue dot under heel centre maximises shock absorption and comfort at point of impact
- Will not deform or degrade with use
- Anti-microbial cover





M-Gel Moisturising Vented Heel Sleeve

| Code | Size |
|--------|----------|
| PT0378 | One Size |



- Gel heel sleeve provides deep hydration to dry, chafed, cracked or damaged skin
- Soft stretch fabric easily rolls on to cover the heel
- New vented fabric design allows for breathability of the skin
- Dermatologist tested



M-Gel Moisturising Booties

| Code | Colour | Size |
|--------|----------------|----------|
| PT0375 | Green | One Size |
| PT0376 | Blue/Turquoise | One Size |



Moisturises and softens the entire plantar surface of the foot as well as heels and toes. The proprietary gel lining releases a combination of skin nourishing hydrators, which help soften and improve the foots appearance.



M-Gel Moisturising Heel Sleeve

| Code | Size | |
|--------|----------|--|
| PT0377 | One Size | |



- Soft stretch fabric easily rolls on to cover the heel and hold anatomically designed gel pad securely in place
- Inside coating of soft mineral oil gel soothes entire heel area on contact
- May be used with creams and heel treatments

placing your order:

tel: +44 (0)151 448 1228

Monday to Friday 8am - 5pm Answering service 5pm - 8am

fax: +44 (0)151 448 1008

email: sales@algeos.com

We aim to respond to all enquiries within 24 hours

web: www.podotech.com

www.algeos.com (24 hour ordering)

mail: A. Algeo Ltd. Sheridan House,

Bridge Industrial Estate, Speke Hall Road, Liverpool, United Kingdom, L24 9BH

delivery:

Same day despatch on all orders placed before 11am (dependent on stock levels). Free delivery on all orders over £500.

payment:

First order customers - payment by credit card or cheque. We accept all major credit cards.











